


## 5-3 High-Precision Measurement of Nuclear Materials – Proposal for the Advanced HKED System –




**Fig.5-8 Schematic view of the advanced HKED system**

The advanced HKED system is composed of an energy recovery linac (ERL), a laser supercavity, and KED/XRF detectors. Mono-energetic X-rays generated by collision of high-energy electrons and laser quanta are used for irradiation of sample solutions. Transmitted and scattered X-rays are detected by the KED and XRF detectors, respectively. The size of the ERL accelerator can be reduced using a 6-loop configuration.

Nuclear material accountancy is an important issue for nuclear safeguards. Generally, destructive analysis with mass spectrometry is used for density determination of nuclear materials. However, non-destructive analysis is desired because of simplicity and rapidity of the analysis, as well as for the reduction of discharged nuclear wastes. Hybrid K-edge densitometry (HKED) is a non-destructive analytical method of nuclear materials dissolved in solution. This technique can be applied for quantification of uranium (U) and plutonium (Pu) concentrations.

In the existing HKED system, bremsstrahlung X-rays from an X-ray tube are used as the incident photon source. The broad energy distribution of the bremsstrahlung X-rays produces background counts due to inelastic scattering that affect the counting precision. However, if mono-energetic X-rays are used, efficient KED and X-ray fluorescence (XRF) measurements are possible through reduction of the background counts. Thus, we have proposed an advanced HKED system with mono-energetic X-rays. Fig.5-8 shows a schematic view of the advanced HKED system.

To evaluate the performance of the advanced HKED



**Fig.5-9 X-ray response spectrum**

A Monte Carlo simulation code was used to investigate the X-ray responses from a sample including U (200 g/l) and Pu (2 g/l). The results for mono-energetic and bremsstrahlung incident X-rays are shown with red and blue lines, respectively. A high-peak intensity of the Pu characteristic X-rays ( $\text{PuK}_{\alpha 1}$ ) can be obtained through efficient measurement of the  $\text{PuK}_{\alpha 1}$  line using mono-energetic X-rays.

system, X-ray responses were investigated using a Monte Carlo simulation code. The result is shown in Fig.5-9. The peak intensity of the Pu characteristic X-rays ( $\text{PuK}_{\alpha 1}$ ) can be enhanced using mono-energetic X-rays. Consequently, both the ratio of the peak to the background count and ratio of the peak to the total count are improved in comparison with those obtained when using bremsstrahlung X-rays. The counting statistics in the region of the U K-absorption edge in the KED spectra are also increased using mono-energetic X-rays.

Analysis of low concentrations of neptunium (Np) is difficult with the existing HKED system because of the high-background counts from the Pu characteristic X-rays. However, with the advanced HKED system, by adjusting the energy of the incident X-rays between the K-absorption edges of Np and Pu, it is possible to measure the Np characteristic X-rays without any influence from the Pu characteristic X-rays.

The present technique is expected to facilitate the non-destructive rapid, and precise density determination of nuclear materials in solutions at nuclear power plants and nuclear reprocessing facilities.

### Reference

Shizuma, T. et al., Proposal for an Advanced Hybrid K-edge/XRF Densitometry (HKED) using a Monochromatic Photon Beam from Laser Compton Scattering, Nuclear Instruments and Methods in Physics Research A, vol.654, issue 1, 2011, p.597-603.